

# Smart CPM Chair

## Product User Manual


### Welcome!

Congratulations on your purchase of our awesome smart CPM chair. To set up your new chair, you will typically need only a 110 V power outlet and a room to put the chair in.


Your package includes the smart chair and an adapter.

### Getting Started

1. plug the chair to the power source
2. Once the chair is powered on, you will see the following screen:



3. Start by touching the “New User” button to get registered.  
The registration window looks like the following:




The image shows a registration window titled "Smart CPM Chair" with a dark blue header bar containing a minimize icon, a maximize icon, and a close icon. The main area has a light orange background. It contains six input fields, each with a label to its left: "Name", "Speed", "Angle", "Starting Angle", "Max Cycles", and "Security Pin". The "Name" field has a small cursor at the beginning. At the bottom, there are two buttons: "Cancel" on the left and "Submit" on the right.

Field Label	Input Type
Name	Text
Speed	Text
Angle	Text
Starting Angle	Text
Max Cycles	Text
Security Pin	Text


Buttons: Cancel, Submit

4. Fill correct information, choose a 4 digits security pin for your profile and touch the button that says submit.

5. Now you should be able to see the following screen:




6. Click on your profile name and you should see the following similar screen:



A screenshot of a mobile application window titled "Smart CPM Chair". The window has a dark blue header bar with the title and standard mobile window controls (minimize, maximize, close). The main content area has a light orange background and displays five configuration parameters in a list-like format. At the bottom, there are three buttons: a back button labeled "<-", a "Use" button, and an "Edit" button.

Name	user
Speed	5
Angle	20
Starting Angle	5
Max cycles	40

7. To start using your preset configurations as shown in the profile, touch the button that says Use. This action will request you put in the security pin of the profile. Once you put in your 4 digits pin and log in, you will see the following screen:




8. Touch the button that best describes your pain level at that moment. Once you do that it will redirect you to the following screen where you can press the button that says "Start Session" to begin using the chair.



9. Once you are finished with your session you can log out of your profile by pressing the button that says "Log Out". This should redirect you back to the welcome screen.

## Adjusting Seat Height

1. Go to Preferences, then touch the button that says “Arrange Seat”. Once you do that you should see the following window:



2. While you are adjusting the seat height, the session will be temporarily paused until you finish seat adjustment. This is normal and your session will resume from where it stopped.